

## KHAYAT | Endodontic Surgical Kit

Part Code: IMKHAYAT1



Dr. Bertrand Khayat, D.D.S., MSD

Dr. Khayat is a leading practitioner and lecturer in Microsurgical Endodontic techniques. His current research interests include new techniques for improved canal preparation. Developed in conjunction with Hu-Friedy, his new hand instrument kit was designed to facilitate these procedures.

PART CODE	DESCRIPTION
IMDIN138S	Large Signature Series® 18 Instrument DIN Cassette, Blue*
10-130-05E	5 Scalpel Handle, European Style, Rounded
SHDPV	Microsurgical Scalpel Handle, Stainless Steel
MM4	Micro Mirror Round 2.5mm
PLGBK6	<b>NEW</b> — Khayat Apical Plugger with Satin Steel® Handle
PFI8ABK6	<b>NEW</b> — Khayat Filling Instrument with Satin Steel Handle
CL856	85 Lucas Surgical Curette, Spoon Shape, Short Satin Steel Handle
PLGVEL56	Velvart Apical Plugger with Satin Steel Handle
TP5070	Perma Sharp® General Tissue Pliers, 15cm/6"
PPR3S6	PR-3 Prichard Periosteal, with Satin Steel Handle
SC13/146	13/14 Columbia Curette, Universal with Satin Steel Handle
SJ34/356	34/35 Jacquette Scaler, Miniature, Posterior, with Satin Steel Handle
CD89/926	89/92 Carver Discoid-Cleoids, with Satin Steel Handle
EXDG16/176	DG 16/17 Endodontic Explorer, with Satin Steel Handle
CRM2	Cawood-Minnesota Retractor
NH5021	Perma Sharp Needle Holders, Castroviejo Curved, 14cm
S14SC	14 LaGrange Double-Curved Super Cut Scissors

All instruments sold in full kit sets are also available for individual purchase. \*Kit includes hinged instrument clips.



## **NEW INSTRUMENTS**

## **KHAYAT** | Endodontic Surgical Instruments

The Dr. Bertrand Khayat Endodontic Surgical Kit includes the new **PFI8ABK6** filling instrument, which has one angled, and one straight spatula end. Clinicians often only use one end on conventional spatulas, but both ends of this instrument are useful.

The straight spatula end is excellent for hard to reach areas in complex surgical procedures.

The **PLGBK6** Plugger, with its 4.5mm and 6.0mm working ends, is ideal for packing retrofilling materials. It was specifically designed for use in deep canal preparation.

Dr. Khayat has given over 200 presentations around the world, and has published the results of his research in many endodontic journals, including the Journal of Endodontics (JOE). Active in many associations, Dr. Khayat maintains his private practice in Paris, France.

## **NEW INSTRUMENTS INCLUDED IN THE KIT**

PART CODE	DESCRIPTION
PLGBK6	Khayat Apical Plugger with Satin Steel Handle
PFI8ABK6	Khayat Filling Instrument with Satin Steel Handle



